

The Queen[®]

Queen of bipolar scissors

Precise and powerful



+44 7483 322776

✉ alan@healthcare25.com

🌐 www.hcare25.ae

📍 H Care 25 Medical & Surgical Trading LLC
16F-A-48, Empire Heights A
Building 16F-A-06
Business Bay, Dubai
United Arab Emirates
Makani No: 2715487271

H CARE25



UAE

The Queen®

Precise. Powerful. Efficient.
Dissection, coagulation, and cutting with a single instrument.

The bipolar scissors The Queen® enable **dissection, coagulation and tissue cutting** in a single step. This reduces instrument exchanges, streamlines the surgical workflow, and shortens procedure times.

One fully insulated blade featuring a **ceramic coating**, together with an insulated connecting screw, ensures precise bipolar current flow for maximum control and safety.

In addition, their **dimensions, design, and handling characteristics** make them an attractive alternative to the **Ethicon® POWERSTAR®** bipolar scissors.



Surface coagulation



Higher efficiency during surgery

- Ceramic-coated technology for controlled energy delivery
- Fewer instrument changes
- Optimized surgical workflow
- Shorter operating times



Pre-dissection coagulation



Efficient bipolar technology

- Reliable and effective hemostasis
- Minimization of thermal stress on surrounding tissue
- Reduced risk of unintended tissue damage
- Bipolar current flows exclusively between the active scissor blades



Simultaneous coagulation and cutting



Precision and Safety

- Insulated scissor blades and connecting screw
- Dissection, coagulation and cutting with a single instrument
- Selected models feature extra-fine tips for microsurgery
- Excellent control when working with delicate tissue structures



Device Compatibility

With four cable options available, bipolar scissors The Queen® ensure a high level of compatibility with most bipolar HF surgical devices.

An overview of available accessories and components can be found at:

www.micromed.com/products/bipolar-scissors/



Precise spot coagulation

Product models

The Queen® bipolar tenotomy scissors, curved, extra-fine tip

143-100-115	115 mm (4 1/2")
-------------	-----------------



The Queen® bipolar tenotomy scissors, curved, extra-fine tip

143-120-175	175 mm (7")
-------------	-------------



The Queen® bipolar Metzenbaum scissors, curved, fine tip

143-310-140	140 mm (5 1/2")
-------------	-----------------

143-320-175	175 mm (7")
-------------	-------------

143-340-230	230 mm (9")
-------------	-------------

143-360-280	280 mm (11")
-------------	--------------



The Queen® bipolar Metzenbaum scissors, curved, standard tip

143-520-175	175 mm (7")
-------------	-------------

143-540-230	230 mm (9")
-------------	-------------

143-560-280	280 mm (11")
-------------	--------------



The Queen® bipolar Mayo scissors, curved, standard tip

143-720-170	170 mm (6 3/4")
-------------	-----------------



Order information

143-100-115	The Queen® bipolar tenotomy scissors, curved, extra-fine tip, 4 1/2"	
		
143-120-175	The Queen® bipolar tenotomy scissors, curved, extra-fine tip, 7"	
		
143-310-140	The Queen® bipolar Metzenbaum scissors, curved, fine tip, 5 1/2"	
		
143-320-175	The Queen® bipolar Metzenbaum scissors, curved, fine tip, 7"	
		
143-340-230	The Queen® bipolar Metzenbaum scissors, curved, fine tip, 9"	
		
143-360-280	The Queen® bipolar Metzenbaum scissors, curved, fine tip, 11"	
		
143-520-175	The Queen® bipolar Metzenbaum scissors, curved, standard tip, 7"	
		
143-540-230	The Queen® bipolar Metzenbaum scissors, curved, standard tip, 9"	
		
143-560-280	The Queen® bipolar Metzenbaum scissors, curved, standard tip, 11"	
		
143-720-170	The Queen® bipolar Mayo scissors, curved, standard tip, 6 3/4"	
		



☎ +44 7483 322776
✉ alan@healthcare25.com
🌐 www.hcare25.ae

📍 H Care 25 Medical & Surgical Trading LLC
16F-A-48, Empire Heights A
Building 16F-A-06
Business Bay, Dubai
United Arab Emirates
Makani No: 2715487271

H CARE25

