



---

# Blepharoplasty

Instrumente für die Blepharoplastik

Instruments for Blepharoplasty

Instrumentos para la blefaroplastia

Instruments pour la blépharoplastie


Strumenti per la blefaroplastica

---


**Instrumente für die Blepharoplastik**  
**Instruments for Blepharoplasty**  
**Instrumentos para la blefaroplastia**  
**Instrumentes pour la blépharoplastie**  
**Strumenti per la blefaroplastica**

Feine Wundhaken mit Isolierungsschicht  
 Fine retractors with insulation layer  
 Separadores finos con aislamiento  
 Ecarteurs fins avec revêtement isolant  
 Divaricatori fini con strato isolante


- |            |                                                                                                                                                                                                       |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| VORTEILE   | <ul style="list-style-type: none"> <li>▶ Verhindert Kontakt zwischen HF-Nadel und Wundhaken</li> <li>▶ Vermeidet Gewebeschädigungen durch punktuelle Hitzeentwicklung</li> </ul>                      |
| ADVANTAGES | <ul style="list-style-type: none"> <li>▶ Prevents contact between HF-needle and retractor</li> <li>▶ Avoids tissue damages caused by punctual heat development</li> </ul>                             |
| VENTAJAS   | <ul style="list-style-type: none"> <li>▶ Evita el contacto entre la aguja de alta frecuencia y el separador</li> <li>▶ Impide del traumatismo de tejido por el desarrollo de calor puntual</li> </ul> |
| AVANTAGES  | <ul style="list-style-type: none"> <li>▶ Empêche tout contact entre l'aiguille HF et l'écarteur</li> <li>▶ Evite le développement ponctuel de chaleur n'endommage les tissus</li> </ul>               |
| VANTAGGI   | <ul style="list-style-type: none"> <li>▶ Impedisce il contatto tra l'ago ad alta frequenza e il divaricatore</li> <li>▶ Evita danni ai tessuti causati dal calore locale</li> </ul>                   |



Fomon



Desmarres



Strandell-Stille



Langenbeck-Green



Die Produkte sind nicht proportional dargestellt.  
 The products are not shown proportionally.  
 Los productos no se muestran proporcionalmente.  
 Les produits ne sont pas représentés de manière proportionnelle.  
 I prodotti non sono rappresentati in modo proporzionale.




**Instrumente für die Blepharoplastik**  
**Instruments for Blepharoplasty**  
**Instrumentos para la blefaroplastia**  
**Instrumentes pour la blépharoplastie**  
**Strumenti per la blefaroplastica**


Blepharopinzette  
Blepharoplasty forceps  
Pinza para blefaroplastia  
Pince de blépharoplastie  
Pinza per blefaroplastica

Messzirkel  
Caliper  
Compás  
Compas  
Compasso

Spezielle Biegung ermöglicht punktuelles Fassen des Gewebes  
Special bending permits punctual tissue grabbing  
Gracias al doblamiento especial el agarre de tejido puntual es posible  
Sa courbure spéciale permet de saisir le tissu de manière ponctuelle  
La curvatura speciale rende possibile la presa precisa del tessuto



Gubisch



Castroviejo

Siebvorschlag  
Set suggestion  
Proposicion de juego  
Proposition pour boîte  
Proposte di set

| Item number | Description                                                             | Quantity |
|-------------|-------------------------------------------------------------------------|----------|
| 04.21.14    | BROPHY (Sullivan) Delicate dissecting scissors · curved · sharp · 14 cm | 1        |
| 11.74.25    | BARRAQUER Micro needle holder · round handle Ø 8 mm · plain · 14 cm     | 1        |
| 20.12.01    | FOMON Retractor 2t ball tips · 16 cm                                    | 1        |
| 20.12.02    | DESMARRES Lid retractor · 14 cm                                         | 1        |
| 20.12.03    | STRANDELL STILLE Retractor · 6 x 17,70 mm · serrated · 19 cm            | 2        |
| 20.12.04    | LANGENBECK-GREEN Retractor · 6 x 25 mm · serrated · 17 cm               | 1        |
| 60.00.40    | CASTROVIEJO Caliper · 20 mm · 8 cm                                      | 1        |
| 60.05.95    | JAEGER Lid plate · 11 cm                                                | 1        |
| 60.33.17    | MOD. BONN Irrigation cannula · 19 Gauge · 3.5 cm                        | 1        |
| 60.34.10    | BOWMAN Lacrimal probe · figure 0000-000 · 13 cm                         | 1        |
| 60.34.31    | BOWMAN (Williams) Lacrimal probe · figure 00-0 · 13 cm                  | 1        |
| 60.34.81    | WILDER Lacrimal dilator · figure 1 · 11 cm                              | 1        |
| • 60.34.82  | WILDER Lacrimal dilator · figure 2 · 11 cm                              | 1        |
| 60.47.05    | CASTROVIEJO Suturing forceps · 0,5 mm · 1 x 2 teeth · 10 cm             | 1        |
| 60.55.45    | GUBISCH Blepharoplasty forceps · 10 cm                                  | 1        |
| 60.73.93    | MOD. BONN Iris scissors · curved · 9 cm                                 | 1        |
| 87.32.12    | Needle electrode · straight · Ø 4 mm                                    | 1        |
| • 87.33.12  | Needle electrode · straight · Ø 2.4 mm                                  | 1        |
| 88.08.01    | Forcep bipolar Non-Stick bayonet · 20 cm                                | 1        |
| 49.50.10    | Medicontainer · 290 x 290 x 100 mm                                      | 1        |
| 49.90.00    | Tray · 190 x 102 x 19 mm · with silicone mat                            | 1        |
| 49.97.90    | Wire basket · 240 x 250 x 40 mm · with silicone mat                     | 1        |



Messbereich  
Measuring range  
Alcance de medición  
Etendue de mesure  
Campo di misura



Instrumente, optional  
Instruments, optional  
Instrumentos, opcional  
Instruments, optionnel  
Strumenti, opzionale



gezahnt  
serrated  
dentada  
dentelée  
dentellato

**Medicon eG**

Gänsäcker 15 · D-78532 Tuttlingen

Tel.: +49 (0) 74 62 / 20 09-0

Fax: +49 (0) 74 62 / 20 09-50

E-Mail: [sales@medicon.de](mailto:sales@medicon.de)Internet: [www.medicon.de](http://www.medicon.de)

Germany



Not all products illustrated in this documentation are available for sale in all countries.  
Not for all declared products in this documentation Medicon is the legal manufacturer.  
Please contact your local distributor for more information concerning availability and regulatory information (e.g. classification) of the mentioned products.

Gedruckt in Deutschland

Modelländerungen vorbehalten

Printed in Germany

Patterns are subject to change

Impreso en Alemania

Salvo modificaciones

Imprimé en Allemagne

Tous droits réservés des changements de modèle

Stampato in R.F.G.

Ci riserviamo la facoltà di cambiamenti nei modelli

Trykkt i Tyskland

453.03.01 © Copyright 06/2025, MEDICON eG, Tuttlingen, Germany